




- General Battery Lifespan is 12 months HOWEVER with the safety precautions within the EquiVibe units and proper battery care, the battery life can be extensively extended.
- Do not store the batteries flat or semi charged for extended periods of time
- Regularly recharge the batteries until fully charged
- Don't leave batteries in direct sunlight or in heat/cold as this will diminish the battery
- Batteries that remain dormant for extended periods of time will reduce capacity and/or stop working.



To replace your battery, take it to your local battery shop, such as Battery World, they are readily available in most places.

Carefully remove the plugs from the old battery and push onto the new battery